FREEVENT®

Neckband





Disclaimer




Atos Medical offers no warranty - neither expressed nor implied - to the purchaser hereunder as to the lifetime of the product delivered, which may vary with individual use and biological conditions. Furthermore, Atos Medical offers no warranty of merchantability or fitness of the product for any particular purpose.

Table 1: Adjust for neck size / Anpassung für die Halsgröße / Aanpassen voor nek grootte / Ajuster pour la taille du cou / regolare per le dimensioni del collo / ajustar para el tamaño del cuello / ajustar para o tamanho do pescoço / Justera för nackstorlek / juster til nakke størrelse / juster for nakke størrelse / Säätää kaulan kokoa / dostosuj rozmiar szyi / prilagodite se veličini vrata

		A-A	B-B	C-C
1651 Freevent Neckband one-piece, small	32-42 cm	24-34 cm	16-26 cm	N/A
1661 Freevent Neckband one-piece, large	45-55 cm	37-47 cm	29-39 cm	21-31 cm

Contents

EN - ENGLISH.....	6
DE - DEUTSCH.....	7
NL - NEDERLANDS.....	8
FR - FRANÇAIS.....	9
IT - ITALIANO	10
ES - ESPAÑOL.....	11
PT - PORTUGUÊS.....	12
SV - SVENSKA	13
DA - DANSK.....	14
NO - NORSK.....	15
FI - SUOMI.....	16
LT - LIETUVIŲ KALBA	17
CS - ČESKY	18
HU - MAGYAR	19
SK - SLOVENČINA.....	20
SL - SLOVENŠČINA	21
PL - POLSKI.....	22
RO - ROMÂNĂ	23
HR - HRVATSKI.....	24
SR - SRPSKI	25
TR - TÜRKÇE.....	26
RU - РУССКИЙ.....	27
AR - العربية.....	28

ENGLISH

Intended use

FREEVENT® Neckband is used for holding a tube or button in a tracheostoma. It is a single use device intended for adult and pediatric patients. Not for neonatal use.

Warnings and precautions

Make sure that the Freevent Neckband is securely fastened to prevent from disconnecting from the tube. To avoid infection please replace the Freevent Neckband every 24 hours. Pressure-induced trauma may appear if the Freevent Neckband is worn too tightly.

Instructions for use

Freevent Neckband two-piece with VELCRO® brand fasteners. Thread both VELCRO® and foam narrow part through the eyelets of the tube neck flange; shortest piece first (1-2). Secure wide VELCRO® brand tab of short section to the long neckband section at proper length and compression. Trim excess material of the long section (3).

Freevent Neckband one-piece with twill ties.

Thread twill ties through the slits in the neckband (1) and trim the neckband to adjust for neck size without cutting the twill ties (Table 1). We recommend to make a knot on the twill before tying it to the tube neck flange (2). Thread one of the twill ties through the neck flange on one side and tie with a double knot (3) to secure the twill to the neck flange. Then do the same on the other side (4). Make sure both twill ties are tied and secured tightly. Trim excess twill material with a straight cut (5) leave some excess to be able to untie.

Disposal

Always follow medical practice and national requirements regarding biohazard when disposing of a used medical device.

Reporting

Please note that any serious incident that has occurred in relation to the device shall be reported to the manufacturer and the national authority of the country in which the user and/or patient resides.

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© Coloplast A/S. All rights reserved.
Atos and the Atos Medical logo are trademarks of Coloplast A/S.
Atos Medical UK Ltd Tottle Road
Atos Medical AB