


# Technical Data



This illustration shows Tracoe Mini Neo tube (REF 350)

## Tracoe Mini for Neonates and Infants

Size	ID	OD	A	C	B	TL	OB	BA
	mm	mm	mm	mm	mm	mm	mm	∅°
2.5	2.5	3.6	6	22	2	30	32	120°
3.0	3.0	4.3	7	23	2	32	34	120°
3.5	3.5	5.0	8	23	3	34	37	120°
4.0	4.0	5.6	9	24	3	36	40	120°

**ID**: inner diameter (clear width) of tube; **OD**: outer diameter of tube;  
**A**: length; **C**: length; **B**: length; **TL**: length along center line from start of neck flange to bottom of tube; **OB**: length along outer bend from start of neck flange to bottom of tube; **BA**: bending angle